

Plattenname:	20260114_VitD	Zeit:	20:42:24	O.D. Obergrenze:	3,500
Anwender:	Fischer	Datum:	14.01.2026	Wellenlängen:	450nm/620nm
Test:	C:\ProgramData\Euroimmun_Analyzer_I\Assays\25-OH Vitamin D.asy (6bd7)				


Kit	E251127AF	261126		
25-OH Vitamin D Biotin + SamBu 1:100		261126		
25-OH Vitamin D S1	0013160223	261126	2,343 O.D.	>1,200 O.D.
25-OH Vitamin D S2	0013170223	261126	1,912 O.D.	
25-OH Vitamin D S3	0013180223	261126	1,505 O.D.	>0,500 O.D.
25-OH Vitamin D S4	0013190223	261126	1,111 O.D.	
25-OH Vitamin D S5	0013200223	261126	0,722 O.D.	>0,100 O.D.
25-OH Vitamin D S6	0013210223	261126	0,397 O.D.	
25-OH Vitamin D PCQ1	0013220223	261126	17,3 ng/ml	8,7-26,0 ng/ml
25-OH Vitamin D PCQ2	0013230223	261126	38,6 ng/ml	27,0-50,2 ng/ml
25-OH Vitamin D Conj	0013540217	261126		
EU Substrate	6006510209			
EU Stop Solution	7006510251			
EU Wash Buffer		261126		

### Validierungskriterien

	S1>S2>S3	2,658>2,156>1,765
	S3>S4>S5	1,765>1,277>0,785
	S5>S6	0,785>0,426
S1	S1>1,200	2,658>1,2
S3	S3>0,500	1,765>0,5
S5	S5>0,100	0,785>0,1

### Quantitative Ergebnisse

Datenmodell 4-Parameter  
A = 2,702                      B = 0,744                      C = 31,45                      D = -0,4102



Plattenname: 20260114\_VitD      Zeit: 20:42:24      O.D. Obergrenze: 3,500  
Anwender: Fischer      Datum: 14.01.2026      Wellenlängen: 450nm/620nm  
Test: C:\ProgramData\Euroimmun\_Analyzer\_I\Assays\25-OH Vitamin D.asy (6bd7)

**Validierungskriterien**

PCQ1      8,7<=PCQ1<=26,0      8,7<=15,4471<=26  
PCQ2      27,0<=PCQ2<=50,2      27<=38,6315<=50,2

**Qualitative Ergebnisse**

If 'Sample>0' Then Result:=\*\*

Standardergebnis := \*

\* No evaluation possible  
\*\*\*\*\* No calculation possible  
MP Manual pipetting  
9,999 OD>3.500  
OD<0.000 PLATE VERIFICATION RECOMMENDED!

**Kombinierter Bericht**

25-OH Vitamin D

Patient ID	O.D. Reader value	Quant. 1	ng/ml value	Well Location	Flag
S1	2,658		0,10	A1	
S2	2,156		3,93	B1	
S3	1,765		10,14	C1	
S4	1,277		25,06	D1	
S5	0,785		59,34	E1	
S6	0,426		>120	F1	

Patient ID	O.D. Reader value	Quant. 1	ng/ml value	Well Location	Flag
PCQ1	1,548		15,45	G1	
PCQ2	1,027		38,63	H1	

Patient ID	O.D. Reader value	Quant. 1	ng/ml value	Well Location	Flag
54001401	0,975		42,30	A2	
24009201	0,756		62,58	B2	
24009001	1,148		31,33	C2	
24009401	0,857		52,10	D2	
75004001	0,865		51,36	E2	
24009101	0,681		72,00	F2	
25005701	0,742		64,22	G2	
64002801	1,106		33,69	H2	
64003001	0,764		61,67	A3	
25006001	0,857		52,10	B3	
44002101	0,912		47,25	C3	
44001801	1,107		33,63	D3	
75003701	1,502		16,81	E3	
34007901	0,903		48,01	F3	
34007401	0,727		66,04	G3	
34007801	0,944		44,66	H3	
45003501	0,841		53,62	A4	
75004101	0,734		65,18	B4	
25005801	0,921		46,50	C4	
55004701	1,117		33,06	D4	
34008001	0,691		70,65	E4	
45003601	1,061		36,42	F4	
34007701	1,294		24,33	G4	
45003401	1,131		32,27	H4	
65004801	0,665		74,23	A5	
75003801	0,762		61,90	B5	
45003801	0,931		45,69	C5	
65004901	0,840		53,72	D5	
35007301	0,710		68,17	E5	
74002301	0,744		63,99	F5	
54002201	1,034		38,16	G5	
74002401	0,874		50,55	H5	
44002001	1,378		20,99	A6	
74002201	1,014		39,51	B6	
64003201	0,922		46,42	C6	
74001901	0,767		61,33	D6	

34008101	1,514	16,44	E6
74002101	1,136	31,99	F6
44001901	0,869	51,00	G6
65004701	1,125	32,60	H6
34008401	1,689	11,81	A7
45003901	1,060	36,48	B7
35007901	1,181	29,59	C7
35007701	1,052	36,99	D7
34008801	0,983	41,71	E7
65005101	0,857	52,10	F7
24010101	1,041	37,70	G7
34008601	0,873	50,64	H7
65005001	0,862	51,64	A8
64003101	0,763	61,78	B8
65005201	0,925	46,18	C8
24010201	0,842	53,52	D8
65006001	0,693	70,39	E8
64003401	0,979	42,00	F8
65005601	0,703	69,07	G8
64003301	0,828	54,89	H8
45004201	0,918	46,75	A9
25006401	0,984	41,64	B9
45004301	0,802	57,53	C9
24010501	1,412	19,77	D9
65005401	0,785	59,34	E9
24010601	0,835	54,20	F9
65006101	0,957	43,65	G9
34009001	0,900	48,26	H9
34009101	0,858	52,01	A10
34009301	0,943	44,74	B10
65005701	1,154	31,01	C10
44002201	0,790	58,80	D10
65005501	0,702	69,20	E10
75004701	1,275	25,15	F10
45004501	0,820	55,69	G10
55004901	0,964	43,12	H10
35008701	0,626	80,02	A11
34009201	0,662	74,66	B11
25006501	0,956	43,73	C11
35008601	0,639	78,03	D11
24010801	0,978	42,08	E11
44002501	0,810	56,70	F11
75004301	0,532	96,51	G11
25006701	1,624	13,39	H11
75004401	0,980	41,93	A12
45004701	1,094	34,40	B12
55005001	0,843	53,43	C12
34007501	0,690	70,79	D12
35008401	0,851	52,67	E12
24011001	0,727	66,04	F12
55005201	0,724	66,41	G12
75004901	1,192	29,04	H12

	Soll		Ist	Temperatur			%OK
	Inkubation	Temperatur	Inkubation	Min.	Mittel	Max.	
1.	35'00	RT-Inkubator	39'12	20,1°C	20,2°C	20,2°C	100
2.	120'00	RT-Inkubator	122'48	20,2°C	20,4°C	20,6°C	100
3.	30'00	RT-Inkubator	29'58	20,4°C	20,4°C	20,6°C	100
4.	15'00	RT-Inkubator	15'01	20,3°C	20,4°C	20,4°C	100